



ColdMark Instructions



STORAGE: For best shelf-life, store ColdMark at room temperature or a cool environment at least 5°C above the response temperature. The indicator is warranted for 2 years from date of sale

ACTIVATION AND APPLICATION: The ColdMark is self-activating, there is no activator tab to pull. Place the indicator on or in your goods using the self-adhesive back. The indicator activates at its stated temperature +/- 1°C

INDICATION: When exposed to temperatures below the response temperatures the bulb irreversibly changes from clear and colorless to a violet color.